

Screenshots

Main Screen

20:07

TDM Insight Vancomycin

Vancomycin TDM

Reset

Enter a case, review the values, then calculate.

Calculation method

Choose the available laboratory samples.

PrePostPre + Post

Sawchuk-Zaske two-point pharmacokinetic estimation using both peak and trough samples.

Load example

Import report

1. Patient information

Patient Name

e.g. Alex Smith

Case ID / MRN

e.g. 123

Age (yrs)

e.g. 52

Gender

MaleFemale

Height (cm)

e.g. 172

Weight (kg)

e.g. 70

Serum Creatinine (μmol/L)

e.g. 85

Unit

μmol/Lmg/dL

Infection Target Profile

Uncomplicated (10-15 mg/L)Severe MRSA (15-20 mg/L)

Calculator

Lab Camera

PK Results

Vancomycin Pre Workflow

20:09

TDM Insight Vancomycin

Vancomycin Pre Workflow

PrePostPre + Post

Calculation based on a pre-dose (trough) concentration.

Load example

Import report

1. Patient information

Patient Name

Puan Fatimah Daud

Case ID / MRN

HRPZ-2026-319

Age (yrs)

74

Gender

MaleFemale

Height (cm)

155.0

Weight (kg)

54.0

Serum Creatinine (μmol/L)

165.0

Unit

μmol/Lmg/dL

Infection Target Profile

Uncomplicated (10-15 mg/L)Severe MRSA (15-20 mg/L)

Loaded: Case 2: Geriatric AKI / Renal Impairment (Pre Only)

2. Dose & timing

Calculator

Lab Camera

PK Results

Vancomycin Post Workflow

20:09

TDM Insight Vancomycin

Calculation method

Choose the available laboratory samples.

Pre

Post

Pre + Post

• Calculation based on a post-dose (peak) concentration and sampling delay.

Load example

Import report

1. Patient information

Patient Name

Chong Wei Lun

Case ID / MRN

HKL-2026-904

Age (yrs)

28

Gender

Male

Female

Height (cm)

180.0

Weight (kg)

82.0

Serum Creatinine (μmol/L)

52.0

Unit

μmol/L

mg/dL

Infection Target Profile

Uncomplicated (10-15 mg/L)

Severe MRSA (15-20 mg/L)

Loaded: Case 3: ICU Polytrauma / High Clearance (Post Only)

3. Dose & timing

Calculator

Lab Camera

PK Results

Vancomycin Pre + Post Workflow

20:10

TDM Insight Vancomycin

Calculation method

Choose the available laboratory samples.

Pre

Post

Pre + Post

• Sawchuk-Zaske two-point pharmacokinetic estimation using both peak and trough samples.

Load example

Import report

1. Patient information

Patient Name

Encik Ahmad Razali

Case ID / MRN

HUKM-2026-881

Age (yrs)

52

Gender

Male

Female

Height (cm)

172.0

Weight (kg)

72.0

Serum Creatinine (μmol/L)

85.0

Unit

μmol/L

mg/dL

Infection Target Profile

Uncomplicated (10-15 mg/L)

Severe MRSA (15-20 mg/L)

Loaded: Case 1: Standard Adult Sepsis (Pre + Post)

3. Dose & timing

Calculator


Lab Camera

PK Results


Calculation Results



Calculation Explanation




Calculation History



What-If Simulation



Camera / OCR



Dark/Light Mode

